



**FENICE S.p.A.**

Via del Lavoro, 1 - 36078 Valdagno (Vicenza) - ITALY  
+39 0445 424888      www.fenice.com

## TECHNICAL DATA SHEET – Leather Care Div.

# ANTI-SQUEAK

### DESCRIPTION

A water-based cream for pigmented and bycast leather.

### HOW TO USE

**Anti-Squeak** is an excellent product designed to remove the squeaking noise made by leather.

**Anti-Squeak** hydrates and revitalises leather to restore flexibility and improve the feel of leather that has become dried out, stiff or squeaky.

Vacuum the leather with a soft-tipped brush. Clean with **Leather Cleaner**. Leave to dry.

Shake the bottle well. Apply **Anti-Squeak** with a clean cloth, wiping with gentle, circular movements. Do not rub aggressively! Leave to dry naturally.

**Pre-test:** Apply the product to a hidden area and only proceed once the required result has been obtained.

PHYSICAL CONSTANTS	Standard Value	Test Method
Appearance	Whitish viscous fluid	MDA 008
Dry content ( % )	13 ± 2	MDA 010/011
pH ( 1:10 )	8.0 ± 1	MDA 004
Density (g/cm <sup>3</sup> )	1.00 ± 0.05	MDA 018

STORAGE	Rev.
Non-dangerous product. Store between +5/40°C	06.2015

### MINIMUM SHELF-LIFE – HANDLING INSTRUCTIONS

The product is guaranteed for its physical and quality characteristics 24 months from the date of production. It has to be kept in its original packaging, well sealed and stored following the conditions in our technical data.

This document is exclusively reserved for professional users.

All information regarding correct use and any hazardous content are available in the safety data sheets.