



FENICE S.p.A.

Via del Lavoro, 1 - 36078 Valdagno (Vicenza) - ITALY
+39 0445 424888 www.fenice.com

TECHNICAL DATA SHEET – Leather Care Div.

MOLD REMOVER

DESCRIPTION

Mold Remover is formulated for the removal of stains caused by the growth of mold on leathers that were exposed to high humidity.

HOW TO USE

Pre-Test: Use protective gloves during application and avoid eye or skin contact. Apply **Mold Remover** to a hidden spot. If there is colour change, loss, or darkening do not use.

Directions: Dampen a soft cloth with **Mold Remover** and gently rub the area to be treated. It might be necessary to repeat the operation after 24 hours. If the mold has been present for a long time, it might be necessary to re-colour the leather. Use **Leather Cleaner** to remove trace residue left by **Mold Remover**.

In case of leathers for car upholstery use our specific **Car Leather Protector**, otherwise use **Leather Protector**.

PHYSICAL CONSTANTS	Standard Value	Test Method
Appearance	Colourless fluid	MDA 008

STORAGE	Rev.
Non-dangerous product. Store between +5/40°C	10.2010

MINIMUM SHELF-LIFE – HANDLING INSTRUCTIONS

The product is guaranteed for its physical and quality characteristics 24 months from the date of production. It has to be kept in its original packaging, well sealed and stored following the conditions in our technical data.

This document is exclusively reserved for professional users.

All information regarding correct use and any hazardous content are available in the safety data sheets.

All our information and suggestions are supplied on the basis of our present knowledge, and are subject to careful tests. However, we cannot accept responsibility for the use of the product, which must be verified and evaluated for its suitability or otherwise by the customer. Declared physical constants are monitored in production and guaranteed by FENICE S.p.A.